PIERENSTEP Feedbacksystem

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
English

Content

Important notes	28	Care and cleaning	40
Safety instructions	28	Contained in the delivery	40
Warnings	29	Classification	40
Some background knowledge	30	Maintenance and technical safety checks	40
Preparations for use	31	Description of the technical safety controls	41
Ascertaining the correct length	31	Measurement technical controls	41
Setting the correct length	31	Guarantee	41
First time feedback system operations	36	Disposal	42
Application	36	Technical data	43
Further functions	37	Dimensions	44
Power saving mode	37	Symbols	45
Charging the battery	38	Labels	46
Technical data and general information	16		

ABBREVIATED USER INSTRUCTIONS

Battery charging

» The PIERENSTEP crutches contain batteries which must be charged so that the feedback system can function. Please use the enclosed charging unit to charge the batteries as and when necessary before using the crutches.

Adjusting the length

» It is important to set the correct length on the crutches. When standing straight on the healthy leg with the arms hanging down by the side, both the handle piece and the wrist joint should be at the same height.

Switching On

» In order to operate the crutches feedback system, first remove either the dummy plug or the battery charger plug. At this point it is important to note that no weight should bear on the crutches.

WALKING WITH THE PIERENSTEP

3 point gait – Step 1


- » Place the crutches and simultaneously the injured leg on the ground while exerting downward pressure on the crutches.

3 point gait – Step 2

- » Continue to exert a constant downward pressure for the whole duration of the step.

3 point gait – Step 3

- » Discontinue downward pressure on the crutches only then when the injured leg is no longer in contact with the ground. The feedback system will indicate if too much or too little pressure has been exerted.

Dear Customer,

We are very pleased that you have decided to purchase a PIERENSTEP crutch system. Please read carefully all the instructions in this manual and take note of the various tips and recommendations. It is important to read these instructions before using for the first time.

Important notes – Please read before using

Purpose of use

The PIERENSTEP crutches are a medical product. During the 3 point gait you are able to control the partial weight bearing level which has been prescribed by your physician or therapist during the course of each step. Incorrect weight bearing levels exerted on the crutches are indicated acoustically by the measurement system and/or by the incidence of vibration. Please do not use these crutches for any other purpose.

Safety instructions

For your own safety, please read carefully and observe the following instructions:

- » If the PIERENSTEP crutches are damaged so that the function or safe use is compromised, please return immediately to the dealer for checking.
- » All maintenance work and checks are to be carried out exclusively by qualified personnel who have received the necessary specialist training in the use of measuring equipment with correct calibration in accordance with the service handbook.
- » Please keep the crutches well away from liquids. The infiltration of liquid into the crutches is to be avoided. If nevertheless liquid has penetrated the crutches, the PIERENSTEP must be thoroughly dried before re-use and must also be subjected to a measurement technical check.
- » Any deviation from these safety measures are the sole responsibility of the user.
- » Please take account of the stated maximum permissible body weight to which the crutches should be subjected.

- » The crutches are to be checked for damage before using.
- » The rubber feet are to be checked regularly for damage and wear and tear.
- » Please have someone accompany you where difficult walking surfaces are involved so that in case of a possible fall help is available.
- » When going up and down stairs use one hand to hold onto the banister and the other hand to hold (in this case a single) crutch. Note: The feedback system of PIERENSTEP is not working while going stairs up or down.
- » In the case of a fall please ensure that the crutches drop to the side so that they do not represent a stumbling block.
- » When walking, keep the crutches close to the body so that they maintain as vertical an alignment as possible.
- » The crutches are only to be used as an aid to walking and should on no account be used for any other purpose (for example to help in the process of standing up, as a tool etc).
- » For reasons of liability it is not permissible to modify the mechanics in any way.
- » Please do not hang articles on the crutches (bags etc.).
- » When the crutches have not been used for longer periods they should be checked by specialist personnel.
- » The crutches should not be used when the ground is covered by ice or snow. Caution! Danger of falls arising from sliding crutches.
- » Be especially careful when crossing loose carpeting, damp tiles or floors.
- » The crutches should not be exposed to the effects of sunrays for longer periods (for example through windows, car windows).
- » Do not lean the crutches on heating or any hot sources.

Warnings

- » Portable and mobile high frequency communication equipment can affect electrical components.
- » Operation in the vicinity (approx. 1m) of short wave or microwave therapy components can cause false measurements leading to incorrect feedback data from the PIERENSTEP measuring system.
- » Each crutch can be loaded with maximum 100 kg.

Some background knowledge

To promote the healing process, for example following an accident or orthopaedic surgery, it can be expedient to subject the relevant extremity as soon as possible to a degree of weight bearing – but not the whole weight. Weight bearing is especially important where the healing of bones is involved because the pressure of the body weight stimulates the growth of bones. Total protection is not therefore desirable and nor is full weight bearing because this would represent too great a strain. It is for this reason that physicians and therapists prefer to prescribe the exercise of partial weight bearing activities. With the help of underarm crutches the patient learns just how much weight should be exerted on the leg or foot which is beneficial for the healing process. At the same time he or she should be mindful of the need to walk properly in order to provide optimum support for the healing process. This is not as easy as it appears at first glance. How much weight is sufficient and how much is too much ? Subjective impressions in this respect are not particularly reliable.

It is for these reasons that the measuring system for the PIERENSTEP crutches has been designed. The system displays precisely how much weight is being exerted on the leg. If the weight level is too great or too small, the user receives an alarm in the form of an acoustic or vibration signal. In this manner the feeling for the correct level of weight bearing is easily acquired and the patient can concentrate on walking correctly.

Important: The proper function of the PIERENSTEP crutches and measuring process is dependent on the clear differentiation between the left-hand and the right-hand crutch. The right-hand crutch is identified by a blue handle piece (2) and a sticker showing „right“. The left-hand crutch is identified by a green handle piece and a sticker showing „left“. Please take care not to interchange the two crutches.

Preparations for use

Ascertaining the correct length

It is very important to ascertain the correct length if the PIERENSTEP crutch system is to function correctly. In the first instance please follow the therapist's instructions. However, you can yourself easily establish the correct crutch length. Stand up straight on the healthy leg and allow one arm to hang down loosely.

If necessary use the other hand for support. Place the crutch adjacent to the hanging arm. The handle piece (2) should now be at the same level as the wrist joint. In case of doubt select a longer setting.

Incorrect length

Correct length

Setting the correct length

In order to set the correct crutch length, please proceed as follows:

1. Loosen the adjusting nut (10).


2. Depress with one hand the interlocking knob (9).

3. Clasp the rubber foot (12) with the other hand.
Caution: Be careful not to clasp the measuring component protective cladding (11).

4. Turn the lower sliding sleeve (7) slightly to the right. The interlocking knob should still be just visible.

5. Now slide the lower crutch sliding sleeve (7) to the desired length. **Caution: Do not push the lower sliding sleeve beyond the lowest catch point 10 (8).**
7. Allow the interlocking knob (9) to engage (to click into place - audibly and visibly).

6.

8. Screw the adjusting nut (10) again tight.

9. Follow the same procedure in respect of the second crutch and the same length. Tip: this is most easily achieved by counting the number of free catch points (8) which are visible above the stop knob (9).

First time feedback system operations

Switch on the feedback system in respect of both crutches. This is effected by removing the blind plug (3) on both crutches or disconnecting the charging unit plug from the charging socket (4) on the front side of the handle piece.

Caution: During the course of switching on the feedback system on no account should the crutches be subjected to weight bearing.

There now follows four acoustic signals and the green LED (5) flashes – at first quickly. Wait for a moment until

the green LED (5) flashing frequency has slowed right down. Only now is the PIERENSTEP feedback system for partial weight bearing control ready for operation; the crutches can again be subjected to weight bearing. In the absence of an acoustic feedback signal the yellow LED (5) flashes and this denotes that your feedback system is not yet ready for operation. Charge the battery with the battery charging unit supplied (please refer to the section „Charging the battery“).

Application

Precise implementation of the three point gait is absolutely essential to ensure the feedback system functions correctly. Your therapist will give you the necessary instructions in this respect.

- » Place the crutches, and at the same time the injured leg, on the ground.
- » Note that at the beginning of the step you must immediately exert the correct amount of force on the

crutches. Caution! At this point most patients already exceed the recommended weight on the injured leg.

- » Continue to exert this same amount of force during the whole duration of the step.
- » Only then when the injured leg is no longer in contact with the ground should the crutches be relieved of weight bearing. The PIERENSTEP feedback system will indicate to you if, during the course of the step, the maximum partial weight bearing on the injured leg has exceeded the tolerance level prescribed by the therapist. Depending on the setting which your therapist has prescribed, feedback takes place in the form of an acoustic signal and/or a vibration alarm in the crutch handle piece (2). In this respect please also take into account the information contained in the patient instruction manual which you have received with your PIERENSTEP crutches.

Further functions

Power saving mode

The crutches switch into power saving mode if they are not used for a period of two minutes. This is indicated by the emission of a double acoustic signal from both crutches. The green LED (5) stops flashing. The crutches can then be re-used by subjecting them to weight bearing for approximately one second.

A short acoustic signal indicates that the crutches are again operational. Please do not hang objects (such as jackets, bags etc.) on the crutches since this would cause interference with the operation of the power saving mode.


The meaning of the various LED displays on the crutches:

LED	MODE
Fast flashing green (2 Hz)	Directly after removing the blind plug or if the blind plug has been removed from only one of the crutches. The crutches are switched on but not ready for operation. The zero compensation takes place.
Slow flashing green (0.5 Hz)	The crutches are ready for feedback system operation and partial weight bearing control.
Flashing yellow – low battery charge level.	The battery must be charged. Neither feedback system nor partial weight bearing control is operational.

Charging the battery

When in regular use, the PIERENSTEP crutches batteries should be charged every two days. Otherwise please charge the battery whenever the yellow LED in the handle piece (2) flashes.

The colour yellow signals that the battery charge level is no longer sufficient to ensure correct crutch operation. Always charge both crutch batteries simultaneously and use the battery charger which has been supplied with the crutches. Other battery chargers can damage the batteries.

Connect the battery charger to a power socket. The operations display on the LED of the charger glows yellow.

Connect the battery charger (13) to the charging socket (4) provided in the crutch handle piece (2). After a few seconds the battery charger LED glows red (14) signalling that charging is in process.

Once the charging process has been completed (after approximately two hours), the battery charger LED glows green and the charger goes into conservation of charge mode.

Meaning of the various LED displays on the battery charger

LED	MODUS
Yellow	Battery not connected
Yellow (glows red after a few seconds)	Battery initialisation and analysis
Red	The battery charging process is in progress.
Green with flashing yellow	Intermediate mode with low charging voltage
Green	Battery is fully charged / charge conservation mode
Red/green flickering	Malfunction

Technical data and general information

Care and cleaning

The PIERENSTEP crutches should be cleaned at regular intervals. A soft, fluff-free cloth perhaps together with a little mild cleaning agent is all that is needed for this purpose. Please do not use any sort of cleaning agent which contains an abrasive medium. When cleaning please ensure that no liquid infiltrates the interior of the crutches. Combustible cleaning agents must have completely evaporated before the crutches can be re-used.

Contained in the delivery

The PIERENSTEP delivery consignment contains the following:

Quantity	REF	Description
1	101528	PIERENSTEP crutches (1 Pair)
2	106346	Battery charger SLG 500-3635
1	451600-0276	User manual
	Version 4	

Classification

The product PIERENSTEP has been designated as Class 1 in accordance with Appendix IX of the Medical Device Directive 93/42/EEC.

Maintenance and technical safety checks

PIERENSTEP crutches may be serviced exclusively by authorised personnel or by the manufacturer. Repairs may be carried out exclusively by trained specialists. Correctly calibrated and appropriate measurement instruments must be used for zero compensation procedures. The responsibility for the carrying out of servicing work rests with the personnel who carry out such work. Unauthorised service and maintenance work compromises both the guarantee and the safe operation of the crutches. The PIERENSTEP crutches must undergo safety and measurement technical checks every twelve months and these checks must be carried out by authorised specialists. To qualify for such authorisation, personnel must have had the necessary education,

possess the necessary knowledge, have gained the required practical experience and be trained with regard to the function and operation of the crutches in order to be able to carry out the necessary checks and controls in a proper manner. Such personnel must be able to carry out their work entirely independently.

Description of the technical safety controls

Sight test:

- » The operating instructions manual to hand.
- » Checking that the equipment is complete.
- » Inscriptions on the labels are legible and complete.
- » Mechanical damage on the crutches (cracks or breakages on the metallic and synthetic parts).
- » Function testing: Low battery charge level, buzzer and vibration alarm.
- » Radio communication with the PC.
- » Carrying out function controls in accordance with the operating instructions and service manual.

Measurement technical controls

- » Crutch calibration at weight bearing levels of 10, 15, 20 and 25 kg.
- » The results of the measurement technical controls and the technical safety controls are to be recorded in a test protocol. The date and the name of the person/company who/which has carried out the controls are to be included in the protocol. In the event that a PIERENSTEP crutch does not function properly or is not operationally safe to use, the crutch must be repaired before re-use.

Guarantee

The PIERENSTEP measuring system is guaranteed for a period of one year from the date of delivery if used in accordance with this manual.

The guarantee does not apply if,

- » The instructions in the operating manual have not been observed,
- » Incorrect operating,
- » Improper use or handling,
- » External intervention by unauthorised persons on the product for reparation purposes,
- » Acts of nature beyond control, for example lightning, strikes etc.
- » Transport damage caused by improper packing when returning the product to the manufacturer/dealer,
- » Lack of maintenance or technical safety controls,
- » Operational normal wear and tear which includes, for example, the rubber feet,
- » Use of non-original spare parts and accessories,

Pierenkemper GmbH is not liable for damages which are caused by deliberate or gross negligence or for injuries to life and limb resulting from usage without the exercise of due care and attention. Statutory guarantee regulations remain unaffected.

Disposal

The manufacturer is committed to take back the PIERENSTEP measuring system and to carry out disposal in a proper manner.

Technical data

Voltage supply	2.4 V – 2 x 1.2 V rechargeable battery
Current take-up	Max. 60 mA in active mode. Approx. 1 mA in stand-by mode
Mass storage	4 Megabyte Flash Memory
Dimensions	150 x 1000 x 50 mm (Width x Height x Depth)
Weight	Approx. 1.2 kg per pair (without accessories)
Maximum load per crutch	100 kg
Ambient conditions	Temperature range 0 °C to 35 °C
Relative humidity	30 % to 90 %

Dimensions

Upper part length: $a = 225 \text{ mm}$

Support angle: $\alpha = 15^\circ$

Lower part length: l

Maximum lower part length: 970 mm

Minimum lower part length: 740 mm

Symbols

Attention: Read accompanying documents, especially user manual!

This symbol follows the year of manufacture for the product.

This symbol follows the address of the manufacturer.

This symbol follows the article number of the product.

This symbol follows the serial number of the product.

This product bears the recycling symbol in compliance with the EU Directive 2002/96/EC. This means that at the end of its working life the product is to be submitted, free of charge, to the communal recycling collection center. Disposal

via the household waste system is not permitted. In this way there is made an active contribution to the protection of the environment.

With the attachment of the CE identification the manufacture declares that the products fulfills all the relevant requirements of the EU directive and has also successfully completed the prescribed conformity evaluation process.


This symbol follows the technical data in respect of the internal battery.

This symbol follows the radio frequency of the electro-magnetic radiation emitted by the product.

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